

Batch Verification

GLOW 70 MG

Independent laboratory documentation

LABORATORY

Janoshik

TESTING PERIOD

24 Jun - 07 Jul 2026

Original reports and raw data reproduced in full.

Verification summary

A consolidated overview of the attached independent laboratory reports. Original Janoshik documents and supporting raw data follow unchanged.

PRODUCT

GLOW 70 mg

BATCH

GL702607

LAB

Janoshik

MEASURED COMPOSITION

Task #195430 / 29 Jun 2026

GHK-CU (GHK CONTENT / COPPER)

68.27 mg (59.83 mg / 8.44 mg)

TB-500 (TB4)

13.20 mg

BPC-157

13.13 mg

ENDOTOXIN

1.886 EU/vial

Task #195431 / 24 Jun 2026

MICROBIAL TESTING

TAMC PASS / TYMC PASS

Task #195433 / 07 Jul 2026

HEAVY METALS

As, Cd, Hg and Pb: not detected

Task #195434 / 29 Jun 2026

DOCUMENT ORDER

1 Composition report and chromatogram

4 Heavy-metals report

2 Endotoxin report

5 Heavy-metals raw data and spectra

3 Microbial testing report

TEST REPORT



E-mail: info@janoshik.com
Web: www.janoshik.com

Task Number #195430

Testing ordered > 14 JUN '26

Sample received > 16 JUN '26

Client OPTide

Sample GLOW

Manufacturer OPTide

Batch GL702607

Sample description >

See picture or pictures attached.

Tests requested >

GLOW (GHK or GHK-Cu/ TB-500/ BPC-157) analysis

Results >



GHK-Cu (GHK content) [Copper Content]	68.27 mg (59.83 mg) [8.44 mg]
TB-500 (TB4)	13.20 mg
BPC-157	13.13 mg

Comments >

Analysis conducted > 29 JUN 2026

Signature >

A handwritten signature in black ink, appearing to be 'Janoshik' or similar, located to the right of the 'Signature >' label.

Verify this test at www.janoshik.com/verify/ with the following unique key

RHLSK9BDZWCJ



SHIMADZU

LabSolutions

Analysis Report


<Sample Information>

Sample Name : 23929-1 GLOW 70mg endo+norm.test+HM
 Sample ID : 015
 Data Filename : 23929-1 GLOW 70mg endo+norm.test+HM_018.lcd
 Method Filename : 251007_MeOH-H2O-TFA_peptidy2_Blend_0.5_+780nm Cu.lcm
 Batch Filename : 20260625_2.lcb
 Vial # : 3-28
 Injection Volume : 0.5 uL
 Date Acquired : 26.06.2026 3:19:50
 Date Processed : 29.06.2026 12:41:23

Sample Type : Unknown
 Acquired by : System Administrator
 Processed by : System Administrator

<Chromatogram>

mAU



<Peak Table>

PDA Ch1 220nm

Peak#	Ret. Time	Area	Conc.	Unit	Name	Area%
1	0.589	6119977	0.00			74.917
2	3.565	1252972	13.13	mg	BPC-157	15.338
3	3.822	796054	13.20	mg	TB-500 (TB4)	9.745
Total		8169003				100.000

TEST REPORT



E-mail: info@janoshik.com
Web: www.janoshik.com

Task Number #195431

Testing ordered > 14 JUN '26

Sample received > 16 JUN '26

Client OPTide

Sample GLOW

Manufacturer OPTide

Batch GL702607

Sample description >

See picture or pictures attached.

Tests requested >

Endotoxin Analysis

Results >



Endotoxin	1.886 EU/Vial.

Comments >

Analysis conducted > 24 JUN 2026

Signature >

A handwritten signature in black ink, appearing to be 'Janoshik' or similar, located next to the 'Signature >' label.

Verify this test at www.janoshik.com/verify/ with the following unique key

S4XKE47HZUXK

TEST REPORT



E-mail: info@janoshik.com
Web: www.janoshik.com

Task Number #195433

Testing ordered > 14 JUN '26

Sample received > 16 JUN '26

Client OPTide

Sample GLOW

Manufacturer OPTide

Batch GL702607

Sample description >

See picture or pictures attached.

Tests requested >

Sterility Testing

Results >



TAMC	PASS
TYMC	PASS

Comments >

Analysis conducted > 07 JUL 2026

Signature >



Verify this test at www.janoshik.com/verify/ with the following unique key

EXAH92RGJPA9

TEST REPORT



E-mail: info@janoshik.com
Web: www.janoshik.com

Task Number #195434

Testing ordered > 14 JUN '26

Sample received > 16 JUN '26

Client OPTide

Sample GLOW

Manufacturer OPTide

Batch GL702607

Sample description >

See picture or pictures attached.

Tests requested >

Heavy Metals Analysis

Results >



As	not detected.
Cd	not detected.
Pb	not detected.

Comments >

Hg: not detected. Although the raw data show low ppm values for heavy metals, these do not necessarily indicate true detection. As shown in the attached graphs, the signals are consistent with electronic noise inherent to the detection system; therefore, the results were concluded as "not detected".

Analysis conducted > 29 JUN 2026

Signature >

A handwritten signature in black ink, appearing to be 'Jano', located to the right of the 'Signature >' label.

Verify this test at www.janoshik.com/verify/ with the following unique key

WWJDBH1ALE3Q

Sample : 23929-1
 Operator: MB
 Comment : Comment
 Group : Hg_Pb_Cd
 Date : 2026-06-26 12:50:58

Measurement Condition

Instrument: EDX-7000 Atmosphere: Air Collimator: 10 (mm) Sample Cup: Polypropylene

Analyte	TG kV	uA	FI	Acq. (keV)	Anal. (keV)	Time (sec)	DT (%)
Cd	Rh 50	640-Auto	#5	0 - 40	22.72-23.52	Live- 51	31
Hg	Rh 50	640-Auto	#5	0 - 40	9.74-10.24	Live- 51	31
Pb	Rh 50	640-Auto	#5	0 - 40	12.38-12.88	Live- 51	31

Qualitative Result

Element: Cd, Hg, Pb

Peak List

Channel	Line	keV	Net Int. (cps/uA)
Cd	CdKa	23.12	0.0005 Q
Hg	HgLa	9.98	0.0001 Q
Pb	PbLb1	12.62	0.0001 Q

Quantitative Result

Analyte	Result	[3-sigma]	Proc.-Calc.	Line	Int. (cps/uA)
Cd	3.8 ppm	L [3.6]	Quan-EC	CdKa	0.0005
Hg	0.2 ppm	L [1.2]	Quan-EC	HgLa	0.0001
Pb	2.8 ppm	L [4.3]	Quan-EC	PbLb1	0.0001

JANOSHIK


Sample : 23929-1
Operator: MB

Group : Hg_Pb_Cd
Comment : Comment


Date : 2026-06-26 12:50:58
01/01




[cps/uA] Cd



[cps/uA] Hg



[cps/uA] Pb




Sample : 23929-1
Operator: MB
Comment : Comment
Group : As_test
Date : 2026-06-26 13:16:00

Measurement Condition

Instrument: EDX-7000		Atmosphere: Air		Collimator: 10 (mm)		Sample Cup: Polypropylene		

Analyte	TG kV	uA	FI	Acq. (keV)	Anal. (keV)	Time (sec)	DT (%)	

As	Rh 50	17-Auto	----	0 - 40	10.30-10.80	Live-	76	31



Peak List

Channel	Line	keV	Net Int. (cps/uA)

As	AsKa	10.54	0.1596 Q

Quantitative Result

Analyte	Result	[3-sigma]	Proc.-Calc.	Line	Int. (cps/uA)
As	0.7 ppm	L [3.2]	Quan-EC	AsKa	0.1596